



WHAT MAKES
VAULTZ UNIQUE

ALL OUR CRYO CHAMBERS
ARE **STATE-OF-THE-ART**
HIGH PERFORMANCE
TECHNOLOGY WITH

REAL ROOM TEMPERATURES

EQUALLY DISTRIBUTED
AND STABLE OVER TIME





REAL WHOLE BODY CRYO

Imagine: in just seconds, you could be on your way to a clear mind. A few minutes of extreme cold stimulation in the Vaultz leaves you feeling refreshed and strengthened, but it also clears your head and frees your mind.

The pleasant feeling of wellbeing lasts for days – you'll be amazed at how the dry, cold air can inspire body, soul and spirit!

KEEP YOUR COOL

www.artofcryo.com



® V1_lux in Alpenresort Schwarz, Austria





MADE IN **GERMANY**

30+ years of
experience

Luxury
design

Cryo
Revolution

COOL, COOLER, VAULTZ

Vaultz® Quality:

All Vaultz are designed and produced at the headquarters of the motherhouse in Germany L&R Kältetechnik GmbH Co. KG

We even do our own programming, have own electric and control cabinet production which allows us to use all potentials of the high-quality components at the optimum!

Of all WBC solutions that origin the motherhouse during the last 27 years, none is off due to technical reasons!

® V1_lux in Alpenresort Schwarz, Austria



**1.5-5
MIN**

Instant effects

Well being ↑
Feeling revitalized ↑
Empowerment ↑

Short term effects

Pain relief ↑
Metabolic efflux ↑
Parasympathetic activity ↑

Soreness ↓
Pro-inflammatory response ↓
Central fatigue ↓
Secondary muscle damage ↓
Muscle tone ↓

Long term effects

Pain ↓
Pro-inflammatory response ↓
Parasympathetic activity ↑
White blood cell activity ↑
Anaerobic threshold ↑

Prof. Dr.sc. med. Winfried Papenfuss
"Power from the Cold", 2016, Germany
Translation by Dr. Julian P. Keogh

@UCryoWellness, Dubai, UAE



LIFESTYLE

Potential effects:



EFFECTIVE **STRESS HANDLING** ●

PROMOTES REGENERATION ●

GENERATES **POSITIVE EMOTIONS** ●

INCREASES PERFORMANCE ●

IMPROVES CENTRAL AND
AUTONOMIC NERVOUS SYSTEM

MENTAL AND PHYSICAL **WELL-BEING** ●

IMPROVES SLEEP BEHAVIOR ●

BOOST **ENERGY AND METABOLISM** ●

BOOST IMMUNE SYSTEM ●

INCREASE **ENDURANCE** ●

IMPROVES HORMONE EXPRESSION ●

INCREASES **LIBIDO** ●


FEEL THIS MOMENT..



Whole-body cryo with real temperatures is effective in giving you a clear mind – and a clear mind gives the soul the freedom to breathe, and helps you return to a true sense of self.


Beauty & Aesthetics. Regular use of whole-body cryo with real temperatures can also have a positive effect on the physical appearance, and can help skin look tighter, brighter, and less wrinkled. Combining the use of real whole-body cryotherapy with various non-invasive and aesthetic treatments offers promising results.

WELLNESS



Possible
BENEFITS OF WBC
with real temperatures

WELLNESS




High-performance cryo-chambers Vaultz are being used in modern fitness studios, health clubs, and other sports facilities that want to offer their customers the ability to train like professionals: sleep management, mental and physical balance as well as pre-cooling (prevention, performance build-up, especially agility and endurance), intermediate cooling for training optimization (shortening the mandatory breaks in advanced training) and post-cooling (improved regeneration) can all be helped during whole body cryo stimulation with real temperatures

SPORT



BENEFITS OF WBC with real temperatures

SPORT



Whole-body cryo has a long history as a therapy in the rehabilitation and medical fields, helping with psychological as well as physical ailments. The feeling of wellbeing that whole body cryo with high performance cryo chambers and real temperatures provide means it continues to enjoy increasing popularity for its health and prevention benefits.

REHABILITATION




EXISTING AND POSSIBLE CENTERS

using **WBC** with real temperatures
for medical, aesthetic, rehabilitation and top class sports

- Medical Stadium (Bordeaux, France)
- Rheumatology Center (Pomorskie Centrum Reumatologiczne, Poland)
- INSEP - the National Institute of Sport, Expertise, and Performance (France)
- Luxury Meditation Retreat Centre (Uluwatu, Bali)
- Cryomedic Center (Châtenoy-le-Royal and Dijon, France)
- Outpatient Rehabilitation center (Rehazentrum Junge, Germany)
- Stade Roland Garros (Paris, France)
- Apex Beauty and Wellness Medical Center (Bangkok and Phuket, Thailand)
- Chenot Weggis, Health Wellness Hotel (Switzerland)
- Alpenresort Schwarz, 5-star Wellness Hotel (Tirol, Austria)
- FFR - French Rugby Federation
- FFT - French Tennis Federation
- FC Bayern Munich, Football Team (Germany)
- Hospitals, for patients, and for employees
- Medical Spas
- Medical Centers
- Beauty and Wellness Clinics
- Headquarters / Corporate Offices
- Lifestyle Clubs
- Physical Therapists
- Recovery Centers
- Wellness Resorts
- Boutique and Luxury Hotels
- Top-class Sport Facilities
- Sport Medicine clinics and centers
- Professional Sport Teams






 x 1 person

 -75°C

Technical characteristics	V1
Cell combination	Single
Width	1 240 mm
Height	2 400 mm
Depth	1 560 mm
Weight - Cryochamber	400 kg
PURE ELECTRICITY DRIVEN	
Machine	Behind or Beside
Cooling system	Two stage Cascade
Cooling capacity	1,5 kW
Cooling time	approx. 1,5 h
Work temperature	-75°C
Weight - machine	530 kg
Interior Design	Elegant
Remote Activation	No
Internet connection	Yes
Control panel	Yes
Music playback	No
Timer setup	Yes





 x 1 person

Stainless stel design
by **Capristo**



Technical characteristics	V1 lux
Cell combination	Single
Width	1 440 mm
Height	2 355 mm
Depth	1 540 mm
Weight - Cryochamber	450 kg
PURE ELECTRICITY DRIVEN	
Machine	Variable
Cooling system	Two stage Cascade
Cooling capacity	2,3 kW
Cooling time	approx. 1,5 h
Work temperature	-85°C
Weight - machine	670 kg
LUX Interior Design by Antonio Capristo	Optional
Remote Activation with mobile phone	Optional
Internet connection	Yes
Control panel	Yes
Music playback	Yes
Timer setup	Yes





2 550 mm

4 200 mm

2 400 mm

2 room version

 -60°C
-110°C

 up to 6 persons

Technical characteristics	V2
Cell combination	Two room version
Width	4 200 mm
Height	2 550 mm
Depth	2 400 mm
Weight - Cryochamber	1 460 kg
PURE ELECTRICITY DRIVEN	
Machine	Separated
Cooling system	3 stage cascade
Cooling capacity	6,9 kW
Cooling time	approx. 4,5 h
Work temperature	-60°C, -110°C
Weight - machine	1 500 kg
Internet connection	Yes
Control panel	Designer control panel made of carbon by capristo.de
Camera in the treatment room	Yes
Music playback in the treatment room	Yes
Intercom in the treatment room	Yes
Timer setup	Yes





2 550 mm

6 000 mm

2 400 mm

3 rooms version

-10°C
-60°C
-110°C

 up to 6 persons

Technical characteristics	V3
Cell combination	Three room version
Width	6 000 mm
Height	2 550 mm
Depth	2 400 mm
Weight - Cryochamber	1 760 kg
PURE ELECTRICITY DRIVEN	
Machine	Separated
Cooling system	3 stage cascade
Cooling capacity	6,9 kW
Cooling time	approx. 4,5 h
Work temperature	-10°C, -60°C, -110°C
Weight - machine	1 710 kg
Internet connection	Yes
Control panel	Designer control panel made of carbon by capristo.de
Music playback in the treatment room	Yes
Intercom in the treatment room	Yes
Camera in the treatment room	Yes
Timer setup	Yes



15 (20)

VAULTZ — HIGH-PERFORMANCE CRYO-CHAMBERS



ART OF CRYO EXCELLENCE

TECHNOLOGY

REAL ultra-low temperatures
Pure cold air
Stable temperatures across the whole body
Equally distributed cold
Development, production, construction, installation and service
Effective use of high-quality components

EFFECTS

Programmable and fully automatic state-of-the-art system
Ready to use in advance
Cold walls - Maximum radiation
Controlled airflow - optimal Convection
Real whole body cryo - including head, neck and shoulders


FINANCE

Low operating cost per application
Marginal cost for an extra treatment is extremely low
Additional cost:
• during product downtimes - minimal • handling - automatic
Heat recovery available
Very low ownership cost per year



UNIQUE APP AND SOFTWARE

For the providers' App works as a CRM platform: our software offers you several exclusive & powerful marketing tools to attract new customers and keep existing ones. You will know who are your customers and what do they want. By creating individual application plans, you bind your customers and deliver unique added value. In addition, the voucher & discount management is available to you.



For the Clients: the software calculates the optimal and individual treatment time and treatment sequence, based on personal physical values, personal thermal perception, treatment aim and automatically create a treatment plan. Clients can see the state of their accounts, can record vital data, can integrate diagnostic data, can easily book or cancel appointments and do much more.

Other features to come: when user arranges the next appointments in V1 lux, she can select favorite music, color, and breath training program.



OTHER PRODUCTS: ART OF CRYO FLOW SYSTEM



OTHER PRODUCTS: ART OF CRYO FLOW SYSTEM

ACTS LIKE OUR "SECOND HEART" FOR A LYMPHATIC SYSTEM

Purpose:

Cardiovascular stimulation

Application:

Treatment and prophylaxis of peripheral circulation problems and cellulitis

During the 30-minute Art of Cryo Flow System application in the so-called space tube, the lower half of the body (from the navel downwards) is acted on with normal and negative pressure. This is not only pleasant but also effective! The treatment method holds an international patent.

Potential benefits: Negative pressure wave treatment is ideal for all who wish for **firmer tissues** on the lower extremities such as thighs, stomach, and bottom. The treatment is also used to **improve weak veins**, spider veins, and blockages. Another target group is a customer with **circulatory disorders** (smokers, diabetics, overweight people) or older customers who are interested in **maintaining good health**. The treatments are also interesting for athletes to **promote faster regeneration** after sports activities or potential for **faster recovery** after liposuction.

Weight, kg	145
Size in cm (l x w x h)	220 x 120 x 110
Mains voltage, V	230
Pressure, mBar	20 - 70
The client:	
- max waist, cm	170
- max weight, kg	180

Physiological response of Intermittent Vacuum Therapy (IVT):

- Increase in arterial blood flow
- Increase in venous return
- Increase in lymph flow (lymphatic drainage)
- Endothelial Stimulation
- Traction



Figure Training



Regeneration
Performance tuning



Health &
Vascular training



Wellness & Vitality

8 PROGRAMS:
4 BASIC / 4 ADVANCED LEVEL

NO CONTRAINDICATIONS





HIGH PERFORMANCE CRYO CHAMBERS

V single Vaultz®

V vario Vaultz®

V unical Vaultz®



art of cryo